



CYNOSURE®

Vivace™

FRACTIONAL RF MICRONEEDLING



## 3D PRECISION - EXPERIENCE THE RESULTS.

The Vivace 3D Precision approach is designed to improve treatment outcomes, utilising insulated needles to achieve optimal results for overall skin rejuvenation, treatment of active acne, acne scars and pore reduction across all skin types.

### INDICATIONS INCLUDE:

- REDUCTION OF FINE LINES AND WRINKLES
- SKIN TIGHTENING
- SCAR REVISION
- REDUCE SPIDER VEINS
- IMPROVE STRETCH MARKS
- REDUCTION OF PORE SIZE
- IMPROVED SKIN TEXTURE
- EVEN SKIN TONE

# Vivace®

## #VIVACESKIN

### BEFORE & AFTER



BEFORE



1 MONTH POST 4 TX

Photos courtesy of Dr Shobhna Singh



BEFORE



1 MONTH POST 4 TX

Photos courtesy of Dr Shobhna Singh



BEFORE



AFTER 3 TX



BEFORE



AFTER 3 TX

Photos courtesy of Dr Shobhna Singh



BEFORE



1 MONTH POST 4 TX

Photos courtesy of Dr Shobhna Singh

# VIVACE - THE ADVANCED MICRONEEDLING RF DEVICE FOR TIGHTENING AND CONTOURING THE FACE, NECK AND BODY.

## Vivace™

Utilising the combination of RF and Microneedling to gently heat the targeted area underneath the epidermis for optimal results, less pain and less downtime than other competitive needling devices.



36 insulated  
gold-tipped  
needles ready  
to deliver results  
at three levels  
of penetration

**The Vivace is FDA cleared and TGA listed, and is positioned as the best-in-class microneedling device - combined with radiofrequency to deliver an incomparable patient experience.**

Setting the Vivace apart from the competition is a number of design elements made to improve patient results and patient comfort.

As opposed to microneedling rollers and pens that can lead to tearing of the skin, pain and variable patient results, the Vivace's needle set contains 36 insulated gold-tipped needles ready to deliver results at three levels of penetration.



The Vivace offers Automatic Visual Confirmation, and red and blue LED light. The robotic depth settings present 31 precise, multi-level adjustments down to a depth of 3.5mm at .1mm increments.

The robotic step motor - only found on the Vivace, creates fast and even puncture rates to the skin. The chevron-pattern for the 70 Watt powered RF distributes heat evenly as it delivers skin tightening across the face, reducing redness, hot spots, burning and discomfort that can occur using other devices.



## VIVACE SCIENCE

### **CIT (Collagen Induction Therapy) by Microneedling multi-level RF**

When applied to the skin, under topical anaesthesia, sterile microneedles are used to create microscopic channels deep into the dermis of the skin that can stay open for up to four to six hours, this stimulates an individual's own body to produce new collagen. The precise delivery of RF energy to 31 varying depths of penetration gives the Vivace exceptional outcomes.

### **WHM (Wound Healing Mechanisms) by multi-level RF energy**

The Vivace microneedles are insulated except the end point (0.4mm). Microneedle electrodes deliver bi-polar energy deep into the dermis without damaging the epidermis. Minimally invasive RF thermal energy creates fractional electrocoagulation in the dermis which facilitates a natural wound healing process. WHM promotes the remodelling of collagen, elastin and wound contraction, which enhances skin laxity.

### **LED (Light Emitting Diode)**

Pulses of light activate skin cells to increase collagen production and moisture levels in the skin, to improve skin tone and texture.



**MINIMAL DAMAGE TO EPIDERMIS – COMFORTABLE TREATMENT  
GREAT RESULTS!**

### **THE DISTINCTION OF VIVACE:**

- Dual frequency (1Mhz & 2Mhz)
- Insulated needle set
- 31 multi-level precise needle depths from .5mm to 3.5mm
- 36 needle set with a small .3mm diameter
- Needles have unlimited one-time use
- High power 70watt Bi-polar RF
- Variable pulse duration of 100ms to 800ms
- Chevron-circuit controlled flow
- Robotic soft soothe glide handpiece for gentle insertion
- Red (stimulates collagen) or Blue (disrupts Acne bacteria) LEDs
- 26.4cm colour touch screen control panel
- Data saving program button
- Automatic Visualisation Confirmation
- Finger switch, foot switch operation and automatic mode

**WE PARTNER  
WITH YOU**



**CUSTOMER**



**ADVANCED  
CLINICAL PARTNER**

- Conducts comprehensive didactic and hands-on clinical training and facilitates staff onboarding
- Provides on-site coaching for patient consultations and assessments
- Maintains customer clinical relations with ongoing support



**BUSINESS  
ACCOUNT MANAGER**

- Helps to facilitate the process and onboarding of product into Practice/Clinic
- Maintains customer relationship on an ongoing basis



**FIELD SERVICE  
ENGINEER**

- Manages on-site installation, testing and monitoring
- Provides on-site technical support with full stock of parts on-hand to minimise machine downtime
- Facilitates software updates on an ongoing basis



**VIVACE  
PRODUCT SPECIALIST**

- Works closely with clinics to support growth of Vivace
- Supports initial onboarding of device with Vivace marketing kit
- Provides customer marketing and promotional tools via Online Client Portal
- Supports PR/Promotional initiatives at local level

Our team of committed experts strive to deliver the optimum experience for our customers, from our best-in-class servicing and clinical training to our dedicated product specialists – we are here to support your business.

**Vivace<sup>®</sup>**  
**#VIVACESKIN**

# PRODUCT SPECIFICATIONS

Model Name	Vivace
Power Input	100/230VAC 50/60Hz
Electric Power Consumption	40VA
Dimensions	380(W) x 1140(H) x340(L)mm
Weight	27Kg
Frequency Output Power Control	1Mhz/2Mhz (High Frequency) PWM
Interface	26.4cm LCD Touch Screen

## CONNECT WITH US

 Cynosure Australia

 Cynosure Australia

 Cynosure Australia

#VIVACESKIN

CYNOSURE®

Contact Cynosure for further information about Vivace

+61 2 9484 4546 | infoaustralia@cynosure.com

Cynosure (Australia) Pty Ltd  
14-16 Suakin St,  
Pymble NSW 2073