



CYNOSURE®

StimSure®

FOR A SCULPTED PHYSIQUE




GIVE YOUR PATIENTS A SCULPTED PHYSIQUE AND PERFECTLY DEFINED CONFIDENCE

Introducing StimSure, state-of-the-art, non-invasive muscle building and toning technology from Cynosure. The new electromagnetic device delivers 24,000 full muscle contractions in just 20-30 minutes for natural-looking, consistent results.

- Builds and tones muscle safely with visible results
- Short treatment time between 20 - 30 minutes, depending on the muscle group being targeted
- 6-8 treatments (2 times a week) for optimum results
- Pre-set programs allow personalised treatment plans
- Marketing support and advice, clinical training and technical support from the experts at Cynosure
- Suitable for most patients

DEFINE PHYSIQUES

Discover how to build and tone muscle with our non-invasive, electromagnetic, muscle stimulation treatment.



StimSure®

SHAPE PHYSIQUES
BUILD YOUR BUSINESS

StimSure®

NON-INVASIVE ELECTROMAGNETIC MUSCLE STIMULATION

- **Longest applicator lifespan of 10,000 treatments (30 minute treatment sessions)**
- **Optimal Tesla per application**, which ensures the required rate of muscle contraction
- **24,000 contractions** during a 20-30 minute treatment session
- **Consistent muscle contraction** resulting in more efficient growth of new protein strands and muscle fibres
- **Pre-set and personalised programs** available for versatile treatments that are tailored to each patient
- Simple to operate with a secure fixing belt
- Single applicator can be used and an applicator arm is available
- Can be operated through light clothing
- **No consumables**
- **Fully delegatable procedure**

THE LATEST IN BODY
CONTOURING TECHNOLOGY
TO GIVE YOU THE COMPLETE
PORTFOLIO AND A COMPETITIVE EDGE.



NO SURGERY, NO DOWNTIME JUST RESULTS

StimSure is the newest addition to the Cynosure portfolio of body contouring products, making it the most complete selection in the industry. Combined with WarmSculpting by SculpSure technology for fat destruction and skin tightening treatments with TempSure, you can now offer your patients a full body shaping package.


Achieving desired and consistent results with heat and cool based technologies alone can be a challenge, but with the Cynosure portfolio you can holistically treat your patients and help them attain their body goals.



HOW DOES STIMSURE WORK?




An electromagnetic field from the StimSure paddles, creates an electric current to engage the muscles.



The electric current stimulates an action potential in the motor neuron cells of the bodies muscles.



The muscle contracts like it would during exercise. The stimulation keeps the muscles contracted tense repeatedly – 'maximal tetanic contraction'.



This consistent muscle contraction results in a more efficient growth of muscle fibres.

**WE PARTNER
WITH YOU**



CUSTOMER



**ADVANCED
CLINICAL PARTNER**

- Conducts comprehensive didactic and hands-on clinical training
- Provides on-site coaching for patient assessments and consultations
- Maintains customer clinical relations with on-site visits and support



**BUSINESS
ACCOUNT MANAGER**

- Helps to facilitate the process and onboarding of product into Practice/Clinic
- Manages customer relationship on an ongoing basis



**FIELD SERVICE
ENGINEER**

- Manages on-site installation, testing and monitoring
- Provides on-site technical support with full stock of parts on-hand to eliminate machine downtime
- Facilitates software updates on an ongoing basis



**STIMSURE
MARKETING**

- Supports initial onboarding of device with StimSure marketing kit
- Provides customer marketing and promotional tools via Online Client Portal
- Supports PR/Promotional initiatives at a local level
- Dedicated social media channel to promote StimSure to patients
- Live clinic finder on our website

Our team of committed experts strive to deliver the optimum experience for our customers, from our best in class servicing and clinical training to our dedicated Marketing team – we are here to support your business.

StimSure®

PRODUCT SPECIFICATIONS

Power Supply	Input 220-240 V, 16 A (Normal power)
Dimensions	542 (L) x 501 (W) x 993 (H) mm
Weight	60kg
Contractions	24,000

CONNECT WITH US

 Cynosure Australia

 Cynosure Australia

 Cynosure Australia



Contact Cynosure for further information about StimSure
+61 2 9484 4546 | infoaustralia@cynosure.com

Cynosure (Australia) Pty Ltd
14-16 Suakin St,
Pymble NSW 2073

© 2020. Cynosure, Inc. All rights reserved. Cynosure and StimSure are registered trademarks of Cynosure. Cynosure owns exclusive rights to photography. Use of photography without written permission of Cynosure is prohibited.

This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. StimSure is only available in selected countries and is not available in the USA. Because Cynosure materials are distributed through websites and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Cynosure representative or write to infoaustralia@cynosure.com StimSure is a non-invasive body contouring system, for a treatment intended for muscular atrophy, with recommended application to strengthen and tighten the abdominal, gluteal and thigh muscles. StimSure is an electrically operated medical device generating an electromagnetic field and is able to contract/stimulate the muscles. Individual results may vary and are not guaranteed. Treatments are not intended for weight-loss or for people who are obese. There are no reports of significant side effects. The post treatment effects that may be experienced include, but are not limited to: Muscular pain, temporary muscle spasm, temporary joint or tendon pain, local erythema or skin redness, premature menstruation and/or lower abdominal cramps. Models for illustrative purposes only and not an actual consumer of the product.