

ULTRAcel Q+BEFORE AND AFTERS












Toning & Large pore (Linear 2.0mm)












Jawline (Linear 4.5mm)

Double chin (Linear 4.5mm)












Double chin (Linear 4.5mm)

Thigh










Axillar

Upper arm

CARTRIDGE OPTIONS





Our team of committed experts strive to deliver the optimum experience for our customers, from our best in class servicing and clinical training to our dedicated marketing team – we are here to support your business.



SLIM HANDPIECE

TO IMPROVE TREATMENT AREA VISIBILITY


Better Accessibility


Improve treatment accuracy with no blind spot

Better Visibility

Up to 65% better visibility of targeted treatment area







COMFORTABLE TREATMENT FOR PATIENTS

DOT

O 100% OF PATIENTS

SAID THE DISCOMFORT

WAS MODERATE TO ZERO


LINEAR

O 100% OF PATIENTS

SAID THE DISCOMFORT

WAS MILD TO ZERO

ADDITIONAL
BENEFITS
REDUCES
FLUSHING AND
REDNESS



PRODUCT SPECIFICATIONS

SYSTEM

Frequency	2MHz, 4MHz, 7MHz
Output Power	0.1 - 3.0J
Dimension	360 x 400 x 980mm
Weight	23.5Kg
CARTRIDGE	
Dot Cartridge	Face: Q 1.5 / QS 2.0 / Q 3.0 / Q 4.5
	Body: Q 6.0 / Q 9.0 / Q 11.0 / Q 13.0
	Spacing: 1 - 2mm (0.1mm/step)
	Length: 5 - 25mm (5mm/step)
Linear Cartridge	Face: Linear 2.0 / Linear 3.0 / Linear 4.5
	Body: Linear 6.0 / Linear 9.0 / Linear 11.0 / Linear 13.0
	Spacing: NA (C/W)
	Length: 25mm
Dimension	87 x 47 x 100mm
Weight	180g
HANDPIECE	
Frequency	2MHz, 4MHz, 7MHz
Output Power	0.1J ~3.0J (0.1J/step)
Dimension	253 x 47 x 105mm
Weight	260g

CONNECT WITH US







Contact Cynosure for further information about ULTRAcel Q+

+61 2 9484 4546 | infoaustralia@cynosure.com

Cynosure (Australia) Pty Ltd 14-16 Suakin St, Pymble NSW 2073